

2X CTAB Buffer Protocol (200ml)

20ml 1M Tris-HCL in dH₂O

56ml 5M NaCl → sol 1

16ml 0.5M EDTA → sol 2

4g CTAB powder

4g PVP

Raise volume upto 200ml

Autoclave

T.E Buffer Protocol (100ml)

1ml of 1M Tris-HCL

0.2ml of 0.5M EDTA (pH-8)

Distilled Water

Autoclave

$$\left(\frac{\text{molarity} \times \text{Mol.Wt} \times \text{Req.Vol}}{1000\text{ml}} \right)$$

$$5\text{M NaCl} = 29.2\text{g} \quad \left(\frac{5 \times 58.4 \times 100}{1000} \right)$$

$$1\text{M Tris-HCL} = 6.057\text{g} \quad \left(\frac{1 \times 157.6 \times 50}{1000} \right)$$

(Tris base= 121.14/ HCL= 36.46)

$$0.5\text{M EDTA} = 18.61\text{g} \quad \left(\frac{0.5 \times 372.24 \times 100}{1000} \right)$$

(121.14+ 36.46= 157.6)