

Gel Electrophoresis Protocol

10ml 50X TAE buffer + 490ml dH₂O = 500ml 1X TAE

100ml 1X TAE buffer in beaker + 1g Agarose powder

Microwave beaker for 1 min with 30 sec intervals until transparent.

Add 6μL EtBr

With stopper enclose the caster and add 100ml sol of 1X TAE+ EtBr

Leave for 20 min to solidify with comb

Remove stoppers and put caster in the apparatus

Fill with 1X TAE buffer

FOR DNA

2μL DNA ladder

4μL Loading Dye + 4μL DNA sample mix on parafilm sheet

Load the sample in wells

Voltage- 85v

Ampere-500mA

Time-35 min