

500ml TAE Buffer (50X)

Add 121g tris base

Add 28.55ml glacial acetic acid

50ml 0.5M EDTA

In 300ml dH₂O (350ml)

Later, make the volume upto 500ml

EDTA buffer

80ml Distilled water

Add 18.61g EDTA

$$\left(\frac{0.5 \times 372.24 \times 100}{1000} \right)$$

1 pellet NaOH

Put on Magnetic stirrer (4-5 min+ 50°C)

Add 10ml dH₂O (now 90ml)

Add 1-2 more pellets of NaOH

Again put on magnetic stirrer (40 min)

Now transparent

Add 10ml dH₂O to make 100ml (Almost 9 and half NaOH pellets)

Now check pH (add 1-2 NaOH pellets to make pH-8)

If pH is more than 8 add 1-2drops of HCL until it becomes 8