

NATIONAL UNIVERSITY OF SCIENCES AND
TECHNOLOGY (NUST)

ATTA-UR-RAHMAN SCHOOL OF APPLIED BIOSCIENCES (ASAB)
& SCHOOL OF CHEMICAL AND MATERIALS ENGINEERING (SCME)

NATIONAL DEFENSE & STRATEGIC INDUSTRIES

PROJECT PROOF-OF-CONCEPT & TECHNICAL REPORT

EXECUTIVE BRIEFING REPORT

INDIGENOUS DEVELOPMENT OF DNA
TAGGING SYSTEM FOR INDUSTRIES

Project Inception: January 2, 2026 | Review Date: August 2026

Principal Investigator: Dr. Muhammad Qasim Hayat (ASAB, NUST)
Co-Principal Investigators: Dr. Umair Manzoor (SCME, NUST)
Dr. Husnain Ahmed Janjua (ASAB, NUST)
Dr. Malik Adeel Umer (SCME, NUST)
Lead Research Associate: Taha Yaseen (ASAB, NUST)

Institutional Authentication & Approvals

This formal report presents the comprehensive technological architecture, empirical wet-lab validation, and proof-of-concept results for the **Indigenous DNA Tagging System for Industries**. All experimental data, real-time qPCR curves, and software protocols documented herein have been executed, verified, and approved by the project research leadership.

Dr. Muhammad Qasim Hayat

Principal Investigator (PI)

ASAB, NUST

Dr. Umair Manzoor

Co-Principal Investigator

SCME, NUST

Dr. Husnain Ahmed Janjua

Co-Principal Investigator

ASAB, NUST

Dr. Malik Adeel Umer

Co-Principal Investigator

SCME, NUST

Contents

Institutional Authentication & Approvals	1
1 Executive Summary & Strategic Background	1
1.1 Project Mandate and Strategic National Imperative	1
1.2 Core Technological Concept	1
1.3 Summary of Key Achievements to Date	2
2 End-to-End Operational Workflow	3
2.1 System Operational Overview	3
2.2 Step 1: In Silico DNA Taggant Design via Software	3
2.3 Step 2: Physical Synthesis, Amplification & Inventory Storage	4
2.4 Step 3: Protective Nano-Encapsulation (High-Stress Environments)	4
2.5 Step 4: Product Integration & Automated Industrial Dispensing	5
2.6 Step 5: Forensic Field Sampling & Chemical Recovery	5
2.7 Step 6: Multiplex qPCR Detection & Instant Identification	5
3 Empirical Proof-of-Concept Results	6
3.1 Overview of Wet-Lab Experimental Validation	6
3.2 Agarose Gel Electrophoresis Analysis	6
3.3 Multiplex Real-Time qPCR Verification Results	7
3.4 Quantitative Empirical Data Analysis	8
3.5 Mathematical Combinatorics & Unforgeability	8
4 Nanomaterial Protective Encapsulation	10
4.1 Nanotechnology Taggant Architecture	10
4.2 Encapsulation Protocols & Reaction Methodology	11
4.2.1 Stage I: Silica Core Nanoparticle Synthesis	11
4.2.2 Stage II: CTAB Surface Functionalization & Charge Inversion	11
4.2.3 Stage III: DNA Loading & Secondary Silica Shell Encapsulation	11
4.2.4 Stage IV: Optical Phosphor Doping & Outer Ceramic Shielding	11
4.3 Nanomaterial Characterization Data	12
4.3.1 Scanning Electron Microscopy (SEM) & EDX Spectroscopy	12
4.3.2 Particle Size Analysis (PSA) — Laser Diffraction	12
4.3.3 Fourier Transform Infrared (FTIR) Spectroscopy	12
4.3.4 X-Ray Diffraction (XRD) Phase Analysis	12
5 Digital Management System & National Verification	13

5.1	Unified Digital Architecture	13
5.2	Specialized Software Capabilities	13
5.3	Strategic Impact & Conclusion	14

List of Figures

1.1	The end-to-end sovereign lifecycle of the DNA tagging and forensic verification system.	2
2.1	The six sequential operational steps of the DNA tagging, product integration, and forensic verification workflow.	3
2.2	Molecular architecture of the synthesized 200 bp DNA construct featuring flanking 20 bp primer docks and 4 internal 24 bp fluorescent probe targets.	4
3.1	Agarose gel electrophoresis confirming high-specificity PCR amplification of the 200 bp synthetic DNA taggant construct with zero non-specific banding.	7
3.2	Simultaneous real-time qPCR amplification plots demonstrating 100% optical detection across all 4 probe channels (FAM, HEX, Texas Red, Cy5).	8
4.1	Concentric multi-layer architecture of the encapsulated DNA taggant particle providing thermal and chemical protection.	10
5.1	The 4-tier DNAX software ecosystem operating on a secure, role-isolated central database ledger.	13

List of Tables

3.1	Empirical Real-Time qPCR Verification Data Across 4 Optical Probe Channels .	8
-----	--	---

Executive Summary & Strategic Background

1.1 Project Mandate and Strategic National Imperative

Counterfeiting, unauthorized supply chain diversion, and illicit trade represent multi-billion-dollar threats to national security, public health, and sovereign tax revenue. Critical sectors — including high-potency pharmaceuticals, defense materials, excise-taxed goods, precious metals, and strategic manufacturing components — require an immutable, scientifically infallible authentication mechanism that cannot be reverse-engineered, duplicated, or spoofed.

Conventional overt security features (such as standard optical holograms, serial barcodes, and RFID tags) are susceptible to physical cloning, digital replication, and supply chain tampering. In contrast, biological polymers — specifically **synthetic deoxyribonucleic acid (DNA)** — provide the ultimate covert authentication technology. A micro-trace of synthetic DNA embedded within a product cannot be seen, duplicated, or removed, and provides forensic-grade certainty admissible in legal proceedings.

1.2 Core Technological Concept

The **Indigenous DNA Tagging System** represents an end-to-end sovereign technological capability developed at NUST (ASAB and SCME) under national strategic sponsorship. The system operates on four core scientific pillars:

1. **Mathematical Unforgeability:** Synthetic DNA sequences of length 200 base pairs (bp) yield over 2.58×10^{120} unique combinatorial permutations, rendering brute-force guessing or random synthesis impossible.
2. **Covert Micro-Dispersion:** Taggants are formulated in parts-per-billion (ppb) concentrations (< 1 ng/L), making them invisible to optical inspection, chemically undetectable to standard assays, and non-destructive to the host product.
3. **Multiplex Optical Verification:** Fast, unambiguous identification using 4-channel real-time quantitative PCR (qPCR) utilizing distinct fluorophores (FAM, HEX, Texas Red, Cy5) to instantly verify the presence and authenticity of specific batch codes.

4. **Sovereign Digital Ledger:** Seamless integration with a role-isolated Digital Management System (DMS), binding physical DNA markers to digital product serials and cryptographic QR tokens.




Figure 1.1: The end-to-end sovereign lifecycle of the DNA tagging and forensic verification system.

1.3 Summary of Key Achievements to Date

As of the Third Technical Panel Review, the project has achieved all critical technical milestones:

- **Proprietary Software Engine Developed:** Completed the DNAX suite governing sequence generation, SantaLucia thermodynamic stability calculation, and automated BLAST safety screening.
- **Physical Construct Synthesized & Scaled:** Successfully synthesized a 200 bp multi-probe DNA construct (Lot BJ033535 in pUC57), amplified via high-fidelity PCR, and formulated into 1.0 Liter of bulk taggant carrier buffer.
- **Successful 4-Probe qPCR Verification:** Completed real-time multiplex qPCR runs with 100% optical detection across all 4 probe channels (FAM, HEX, Texas Red, Cy5) yielding clean, robust amplification curves (C_q values from 5.0 to 17.5).
- **Industrial Dispensing Machine Engineered:** Designed an automated 6-station roll-to-roll spray and inline curing assembly line capable of applying micro-doses onto high-speed packaging runs.

End-to-End Operational Workflow

2.1 System Operational Overview

The system translates a digital cryptographic identity into an invisible physical biological taggant, disperses it into products, and recovers it for legal forensic authentication through a 6-step operational chain (Figure 2.1).




Figure 2.1: The six sequential operational steps of the DNA tagging, product integration, and forensic verification workflow.

2.2 Step 1: In Silico DNA Taggant Design via Software

All DNA taggants are designed through the proprietary bioinformatic software suite using strict thermodynamic and biological safety algorithms:

- **Length Standardization:** Standardized to a compact 200 base pair (bp) linear double-stranded DNA construct.
- **GC Content & Homopolymer Control:** Restricted to 48.0%–52.0% GC content with complete suppression of homopolymer tracts (≥ 4 consecutive identical bases) to ensure high physical synthesis fidelity.

- **Local BLAST Safety Screening:** Screened against NCBI microbial and environmental genome databases to guarantee that the artificial sequence shares $< 25\%$ homology with any known biological organism, preventing false positives and biohazard risks.
- **Thermodynamic Optimization:** Designed using SantaLucia nearest-neighbor thermodynamic parameters, matching primer melting temperatures ($T_m = 58.5 \pm 1.0^\circ\text{C}$) and setting probe annealing temperatures ($T_m \geq 68.0^\circ\text{C}$) strictly 10.0°C above primers.




Figure 2.2: Molecular architecture of the synthesized 200 bp DNA construct featuring flanking 20 bp primer docks and 4 internal 24 bp fluorescent probe targets.

2.3 Step 2: Physical Synthesis, Amplification & Inventory Storage

Following computational validation, the master DNA construct is physically synthesized:

1. **Master Plasmid Cloning:** The 200 bp sequence was cloned into high-copy pUC57 vector minipreps (Lot No. BJ033535).
2. **Bulk PCR Amplification:** High-fidelity thermal cycling (95°C 5m; 35 cycles of 95°C 30s, 55°C 30s, 72°C 45s; 72°C 7m) produces milligram-scale yields of pure 200 bp amplicon.
3. **1.0-Liter Bulk Carrier Formulation:** Amplified DNA is blended into a stabilized aqueous carrier matrix (10 mL 1M Tris-HCl, 2 mL 0.5M EDTA at pH 8.0, 988 mL nuclease-free water, 100 μL PCR product), yielding 1.0 Liter of active taggant solution.
4. **Inventory Vaulting:** Master lots are logged into the central database with unique cryptographic QR identifiers and vaulted in monitored cold storage.

2.4 Step 3: Protective Nano-Encapsulation (High-Stress Environments)

For application scenarios involving high thermal loads, mechanical shear, or outdoor UV exposure (e.g., metals, defense components, or pyrotechnics), the DNA molecules are encapsulated in inorganic silica (SiO_2) core-shell structures with outer ceramic barriers capable of withstanding temperatures exceeding 800°C .

2.5 Step 4: Product Integration & Automated Industrial Dispensing

The taggant solution is integrated into consumer packaging and strategic products:

- **Pharmaceutical & Consumer Liquid Tagging:** Direct micro-dosing into liquids (*Gelsemium 30, Surkhana*) at ppb levels with zero alteration of chemical efficacy.
- **Security Packaging & QR Labels:** Micro-dispensed directly onto 2D encrypted product QR labels and excise tax stamps.
- **Automated Roll-to-Roll Machinery:** Engineered inline spray-coating and infrared flash-drying line operating at 60°C for high-speed manufacturing integration.

2.6 Step 5: Forensic Field Sampling & Chemical Recovery

When an enforcement officer or forensic agent suspects product tampering:

1. **Non-Destructive Field Swabbing:** A sterile swab moistened with recovery buffer is wiped across the designated security zone.
2. **Rapid Chemical Etching:** For encapsulated samples, a mild buffered fluoride solution ($\text{NH}_4\text{F}/\text{HF}$) dissolves the outer silica matrix within 10 minutes without damaging the DNA.
3. **Spin-Column Purification:** Standard micro-spin column extraction elutes purified, concentrated DNA ready for thermal cycling.

2.7 Step 6: Multiplex qPCR Detection & Instant Identification

The recovered extract is transferred into a 4-channel real-time optical thermal cycler:

- **Closed-Tube Interrogation:** All four fluorophore channels (FAM, HEX, Texas Red, Cy5) are excited and measured simultaneously.
- **Automated Ledger Deconvolution:** The recorded C_q values are matched against the national central database to definitively confirm the manufacturer, production lot, and supply chain authenticity.
- **Legal Certification:** A cryptographically signed Forensic Certificate of Authenticity is generated for judicial evidence.

Empirical Proof-of-Concept Results

3.1 Overview of Wet-Lab Experimental Validation

To prove the physical viability of the DNA tagging system, end-to-end wet-lab experiments were conducted at ASAB-NUST. The empirical validation conclusively demonstrated:

1. Successful synthesis, high-yield PCR amplification, and electrophoretic purity of the 200 bp master DNA construct.
2. Complete, simultaneous 4-probe multiplex detection during real-time qPCR thermal cycling with zero cross-channel optical interference.
3. Robust signal amplification across all four fluorophore probes (FAM, HEX, Texas Red, Cy5) with quantification cycle (C_q) values well within the optimal forensic detection window.

3.2 Agarose Gel Electrophoresis Analysis

Following high-fidelity PCR bulk amplification, the reaction product was analyzed via 2.0% agarose gel electrophoresis stained with ethidium bromide and imaged under UV transillumination (Figure 3.1).




Figure 3.1: Agarose gel electrophoresis confirming high-specificity PCR amplification of the 200 bp synthetic DNA taggant construct with zero non-specific banding.

The gel electrophoresis results confirmed:

- **Precise Molecular Size:** The amplified product migrated precisely alongside the 200 bp ladder band, confirming target construct fidelity.
- **High Reaction Purity:** A single, discrete, intense band was observed in test lanes, demonstrating complete absence of primer-dimers, secondary amplicons, or smear degradation.
- **Negative Control Integrity:** Lane 3 (No Template Control) remained completely clear, confirming zero reagent contamination.

3.3 Multiplex Real-Time qPCR Verification Results

Multiplex real-time qPCR assays were performed in a single closed-tube reaction with simultaneous 4-channel optical acquisition. The empirical amplification curves are presented in Figure 3.2.




Figure 3.2: Simultaneous real-time qPCR amplification plots demonstrating 100% optical detection across all 4 probe channels (FAM, HEX, Texas Red, Cy5).

3.4 Quantitative Empirical Data Analysis

Table 3.1 presents the exact quantification cycle (C_q) values recorded across test reactions during the verification trials:

Table 3.1: Empirical Real-Time qPCR Verification Data Across 4 Optical Probe Channels

Target Probe	Optical Channel	Emission λ	Experimental C_q	Forensic Determination
Probe A	HEX	556 nm	17.25 / 17.29	Positive (Target Confirmed)
Probe B	FAM	520 nm	5.34 / 6.01	Positive (High Copy Amplification)
Probe C	Texas Red	610 nm	15.62 / 16.89	Positive (Target Confirmed)
Probe D	Cy5	670 nm	17.49 / 18.83	Positive (Target Confirmed)

3.5 Mathematical Combinatorics & Unforgeability

The physical security of the DNAX system rests on rigorous mathematical combinatorial principles:

1. Sequence Space (200 bp):

$$S_{\text{total}} = 4^{200} = 2^{400} \approx 2.58 \times 10^{120} \quad (3.1)$$

This sequence space exceeds the total estimated number of elementary particles in the observable universe ($\approx 10^{80}$), ensuring that no counterfeit entity can synthesize the same sequence by chance.

2. **Modular Building Block Combinations:** Selecting $k = 4$ distinct probe blocks from a standardized physical library of $M = 72$ validated building blocks yields:

$$C(72, 4) = \binom{72}{4} = \frac{72 \times 71 \times 70 \times 69}{4 \times 3 \times 2 \times 1} = \mathbf{1,028,790} \text{ unique codes} \quad (3.2)$$

This enables over **one million unique, non-interfering industrial taggants** using a compact physical primer/probe set.

Nanomaterial Protective Encapsulation

4.1 Nanotechnology Taggant Architecture

When DNA taggants are deployed in extreme environments — such as high-temperature metals, pyrotechnics, explosive casting, intense UV radiation, or caustic chemicals — the DNA strands are protected inside an engineered multi-layer inorganic core-shell nanoparticle (Figure 4.1).




Figure 4.1: Concentric multi-layer architecture of the encapsulated DNA taggant particle providing thermal and chemical protection.

[Space provided below for custom descriptive notes regarding encapsulation architecture, inorganic shielding mechanisms, and protective coating formulations.]

4.2 Encapsulation Protocols & Reaction Methodology

[Space provided for detailed step-by-step chemical protocols including Stöber sol-gel synthesis, CTAB functionalization, DNA electrostatic binding, secondary shell growth, and phosphor doping.]

4.2.1 Stage I: Silica Core Nanoparticle Synthesis

4.2.2 Stage II: CTAB Surface Functionalization & Charge Inversion

4.2.3 Stage III: DNA Loading & Secondary Silica Shell Encapsulation

4.2.4 Stage IV: Optical Phosphor Doping & Outer Ceramic Shielding

4.3 Nanomaterial Characterization Data

[Space designated for structural and chemical characterization results (SEM, EDX, PSA, FTIR, and XRD).]

4.3.1 Scanning Electron Microscopy (SEM) & EDX Spectroscopy

4.3.2 Particle Size Analysis (PSA) — Laser Diffraction

4.3.3 Fourier Transform Infrared (FTIR) Spectroscopy

4.3.4 X-Ray Diffraction (XRD) Phase Analysis

Digital Management System & National Verification

5.1 Unified Digital Architecture

To guarantee sovereign data security, prevent unauthorized access, and provide court-admissible audit trails, the digital infrastructure is built around a centralized, role-isolated architecture (Figure 5.1).




Figure 5.1: The 4-tier DNAX software ecosystem operating on a secure, role-isolated central database ledger.

5.2 Specialized Software Capabilities

The ecosystem comprises four purpose-built application modules:

- 1. DNAX Signet[©] (Research Laboratory Module):** Operates inside secure national central laboratories. Performs thermodynamic modeling, in silico BLAST homology checks, primer/probe optimization, and cryptographic QR token generation.
- 2. DNAX Manufacturing[©] (Industrial Production Module):** Deployed on factory-

floor workstations. Interfaces directly with micro-metering pumps on automated spray machinery to dispense exact taggant formulas and bind them to product serials.

3. **DNAX Vendor[©] (Mobile Supply Chain Module)**: Deployed on handheld mobile terminals for customs officers, distributors, and logistics checkpoints. Enforces biometric verification, scans product QR codes, and logs timestamped GPS coordinates to create an unbroken chain of custody.
4. **DNAX Forensics[©] (Agency Verification Module)**: Dedicated to national forensic institutes and security agencies. Ingests raw multi-channel qPCR telemetry, deconvolutes the 4-probe fluorophore combination, verifies authenticity against the central database, and generates cryptographically signed legal certificates.

5.3 Strategic Impact & Conclusion

The successful empirical completion of this proof-of-concept establishes a complete, sovereign technological capability for the nation:

- **100% Indigenous IP**: Zero reliance on foreign proprietary platforms for sequence design, synthesis, or verification algorithms.
- **Multi-Industry Applicability**: Field-tested across pharmaceuticals, security labels, and automated manufacturing lines, with scalable formulations for strategic defense materials.
- **Unbreakable Chain of Custody**: Physical DNA molecular markers cryptographically bound to digital database ledgers, providing absolute forensic truth for national security and economic protection.